



Getein
Biotech, Inc.

Stock Code
603387

FIA8600

Quantitative Immunoassay Analyzer

**OPTIMIZED SOLUTION TO ACUTE
AND CHRONIC DISEASES**

Specifications

Methodology	Quantitative Immunoassay
Sample	Whole Blood/ Serum/Plasma/ Uri
Language	English/Chinese
Display	7 inch touch screen
Weight	1.5 KG
Printer	Internal Thermal Printer
Data Storage	10,000 data
Dimensions	220 mm×230 mm×98 mm (D×W×H)
Certificate	CE, ISO
Application	Hospital, Clinic, Laboratory
Communication port	USB / LAN



FIA8600 Test Items

Cat.#	Test Item	Disease	Measuring Range	Sample	Cut-off Value	Reaction Time	Qualification
CG3001	cTnI	Myocardial infarction	0.50-50.00ng/mL	S/P/WB	0.50ng/mL	15min	CE NMPA
CG3002	NT-proBNP	Heart failure	100-35000pg/mL	S/P/WB	300pg/mL	15min	CE NMPA
CG3003	hs-CRP+CRP	Cardiovascular inflammatory diseases/ inflammatory disorders	0.5-200.0mg/L	S/P/WB/ Fingertip Blood	10.0mg/L 3.0mg/L	90s	CE NMPA
CG3004	NT-proBNP /cTnI	Heart failure/ acute coronary syndrome	100-12000pg/mL 0.50-50.00ng/mL	S/P/WB	300pg/mL 0.50ng/mL	18min	CE NMPA
CG3005	CK-MB /cTnI /Myo	Myocardial injury	2.50-80.00ng/mL 0.50-50.00ng/mL 30.00-1000.00ng/mL	S/P/WB	5.00ng/mL 0.50ng/mL 70.0ng/mL	15min	CE CFDA
CG3006	D-Dimer	Venous thromboembolism/ pulmonary embolism	0.10-10.00mg/L	S/P/WB	0.50mg/L	7min	CE NMPA
CG3007	PCT	Sepsis/ septic shock	0.10-50.00ng/mL	S/P/WB	0.10ng/mL	15min	CE NMPA
CG3008	CysC	Early diagnosis of kidney disease/ detection of kidney damage for surgery patients	0.50-10.00mg/L	S/P/WB	0.51-1.09 mg/L	3min	CE NMPA
CG3009	mAlb	Early diagnosis and evaluation of diabetic nephropathy	10.0-200.0mg/L	Urine	20.0mg/L	3min	CE NMPA
CG3010	NGAL	The best indicator of early renal injury	50.0-5000.0ng/mL	S/Urine	Serum:200.0ng/mL Urine:100.0ng/mL	3min	CE NMPA
CG3011	β_2 -MG	Kidney damage for diabetic & hypertensive patients	0.50-20.00mg/L	S/P/WB	0.80-3.00 mg/L	3min	CE NMPA
CG3012	CK-MB /cTnI	Myocardial injury	2.50-80.00ng/mL 0.50-50.00ng/mL	S/P/WB	5.00ng/mL 0.50ng/mL	15min	CE
CG3013	HCG+ β	Pregnancy early test	5.0-100000.0mIU/mL	S/P	5.10mIU/mL	10min	CE NMPA
CG3017	HbA1c	Diabetes mellitus	2.00%-14.00%	WB	3.80%-5.80%	3min	CE NMPA NGSP IFCC
CG3018	CK-MB	Myocardial injury	2.50-80.00ng/mL	S/P/WB	5.00ng/mL	15min	CE
CG3024	TSH	Thyroid diseases	0.10-50.00 μ U/mL	S/P	0.27-4.20 μ U/mL	15min	CE
CG3047	<i>H. Pylori</i>	<i>H. Pylori</i> infection	5-200ng/mL	Stool	100 μ L	10min	CE
CG3042	FOB	Gastrointestinal diseases	100-2000ng/mL	Stool	100 μ L	10min	CE

GP Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508/68568594
Fax: +86-25-68568500
E-mail: sales@getein.com.cn; overseas@getein.com.cn
Web: www.getein.com

ISO 9013485 FSC CE NMPA NGSP IFCC IVD

